

Plant Virus

Real-time RT-PCR

Target virus

BPeMV, CFMMV, CGMMV, CMMoV, HLSV, KGMMV, MarMV, ORSV, PMMoV, RMV, SFBV, TMGMV, TMV, ToMV, ZGMMV, WMV, CMV, ZYMV, MNSV, PRSV

Methodology

one-step Real-time RT-PCR

Short & Simple Protocol

After RNA Ext.

Step 1 : Mix extracted RNA with the Kit's reagents

Step 2 : One-step Real-time RT-PCR run

Step 3 : Analysis

2.5 hrs


Kit Components

- Virus specific primer/probe mixture
- 2X RT-PCR Reaction Buffer
- RT enzyme mix
- Control

Storage

-20°C 보관

Exp.



Ordering Information

PowerChek™ Plant Virus Real-time PCR Kit series	Cat #	rxn	PowerChek™ Plant Virus Real-time PCR Kit series	Cat #	rxn
PowerChek™ PWN Real-time PCR Kit (소나무재선충)	R0901	50	PowerChek™ RMV Real-time PCR Kit	R0913	50
PowerChek™ Tobamovirus Real-time PCR Kit	R0902	50	PowerChek™ SFBV Real-time PCR Kit	R0914	50
PowerChek™ RSV Primer/Probe	R0903	100	PowerChek™ TMGMV Real-time PCR Kit	R0915	50
PowerChek™ BPeMV Real-time PCR Kit	R0904	50	PowerChek™ TMV Real-time PCR Kit	R0916	50
PowerChek™ CFMMV Real-time PCR Kit	R0905	50	PowerChek™ ToMV Real-time PCR Kit	R0917	50
PowerChek™ CGMMV Real-time PCR Kit	R0906	50	PowerChek™ ZGMMV Real-time PCR Kit	R0918	50
PowerChek™ CMMoV Real-time PCR Kit	R0907	50	PowerChek™ WMV Real-time PCR Kit	R0919	50
PowerChek™ HLSV Real-time PCR Kit	R0908	50	PowerChek™ CMV Real-time PCR Kit	R0920	50
PowerChek™ KGMMV Real-time PCR Kit	R0909	50	PowerChek™ ZYMV Real-time PCR Kit	R0921	50
PowerChek™ MarMV Real-time PCR Kit	R0910	50	PowerChek™ MNSV Real-time PCR Kit	R0922	50
PowerChek™ ORSV Real-time PCR Kit	R0911	50	PowerChek™ PRSV Real-time PCR Kit	R0923	50
PowerChek™ PMMoV Real-time PCR Kit	R0912	50			